

# Phonak

# RemoteControl

## Technical Data

Dimensions (L x W x H)	84 x 32 x 11.5 mm (3.3 x 1.3 x 0.5")
Weight	26.4 g (9.3 oz), battery included
Color	Champagne/black
Operating conditions	0° to +45° Celsius (32° to 113° Fahrenheit)
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)
Humidity during storage	0% to 70% (non condensing)
Humidity during operation	0% to 95% (non condensing)
Included parts	Lanyard, CR 2032 battery
Ingress protection	IP52

## Power supply

Battery	CR 2032
Operating time	12 month at average use

## Radio information

Antenna type	Monopole antenna
Operation frequency	2.4 GHz – 2.48 GHz
Modulation	GFSK
Radiated power	< 2.8mW



## Standards applied

---

### Emission standards

EN 60601-1-2:2015  
IEC 60601-1-2:2014  
EN 55011:2009+A1  
CISPR11:2009/AMD1:2010  
CISPR22:1997  
CISPR32:2012  
ISO 7637-2:2011  
CISPR25:2016  
EN 55025:2017

---

### Immunity standards

EN 60601-1-2:2015  
IEC 60601-1-2:2014  
EN 61000-4-2:2009  
IEC 61000-4-2:2008  
EN 61000-4-3:2006+A1+A2  
IEC 61000-4-3:2006+A1+A2  
EN 61000-4-4:2012  
IEC 61000-4-4:2012  
EN 61000-4-5:2014  
IEC 61000-4-5:2014  
EN 61000-4-6:2014  
IEC 61000-4-6:2013  
EN 61000-4-8:2010  
IEC 61000-4-8:2009  
EN 61000-4-11:2004  
IEC 61000-4-11:2004  
IEC 60601-1 (§ 4.10.2):2005  
ISO 7637-2:2011

---

### Legal manufacturer



Sonova AG  
Laubisrütistrasse 28  
CH-8712 Stäfa  
Switzerland

---